







スイングロータ Swing Rotor T4SP

最高回転速度: 3,600rpm
 Max. Speed
 最高遠心加速度: 2,028×g
 Max. R.C.F.
 全容量: 90ml = 15ml × 8
 Nominal Rotor Capacity



適用バケット・チューブ・ボトル一覧表 Summary of Buckets, Tubes, Bottles

バケット Bucket 品名 ・ Part No. ・ 寸法 ・ 材質 Name ・ Part No ・ Size ・ Material	チューブ Tube/ボトル Bottle 品名 ・ Part No. ・ 寸法 ・ 実容量 ・ 材質 Name ・ Part No ・ Size ・ Actual capacity ・ Material	最高回転速度 Max. Speed (遠心加速度) (Max. R.C.F) (Rmax)
 標準付属品 Standard accessory 材質 Material: ポリプロピレン Polypropylene	 15ml プラスチックチューブ Plastic Tube (市販品) (Marketed) (D)17X(L)106 ※2	※1 3,600rpm (2,028xg) (Rmax:14.0cm)
 標準付属品 Standard accessory 材質 Material: ポリプロピレン Polypropylene	 7ml プラスチックチューブ Plastic Tube (市販品) (Marketed) (D)13.6X(L)106 ※2	
単位 Unit (mm) 		

※1: 市販品のチューブ、ボトル等は、メーカー指定の許容回転速度または許容遠心加速度以下で使用してください。許容値を超えて使用するとチューブ、ボトル等が破損する場合があります。

※1: When using tubes, bottles, etc. on the market, perform operation under the allowable speed or the allowable RCF specified by the manufacturer and do not use any tube, bottle, etc. whose allowable speed or RCF is not specified. Otherwise the tubes, bottles, etc. may be broken during operation.

※2: 市販プラスチックチューブの実容量は、各製品の仕様により異なりますので、各メーカーのカタログをご覧ください。

※2: Actual capacities of plastic tubes on the market may vary depending on the specifications. Refer to the manufacture's catalog.